


MS Spectrum



Sample Information

Dissolution method	:15% ACN+85% H2O	Interface	:ESI	Prerod Bias	:+4.5kv
Date Acquired	:2025/07/18 16:33:10	Nebulizing Gas Flow	:1.50L/min	Detector	:-0.2kv
Injection Volume	:0.2ul	CDL Temp	:250C	T.Flow	:0.2ml/min
Block Temp	:200	CDL Volt	:0v	B.conc	:50% H2O/50% MEOH

Name	:N/A
Sequence	:G(N15)PA(N15)G(N15)PQG(N15)PR
Theoretical	:839.91
Observed	:839.80